

Svetlana Gramatikova¹, Andrew Carson¹, Bradley Patay¹, Tim Stenzel¹, and Jeffrey E. Miller¹
¹Invivoscribe, San Diego, USA


Abstract

Introduction: Molecular abnormalities in multiple genes cause acute myeloid leukemia (AML), a clinically and genetically heterogeneous disease. A recent study has shown that up to one-third of AML patients would have had their treatment altered using a knowledge bank of matched genomic – clinical data, as compared to the treatment they received using current practice recommendations¹. With a frequency of ~30%, nucleophosmin (*NPM1*)-mutated AML is the largest genomic subgroup in this disease^{2,3}. The key purpose of this study is characterization of the genetic heterogeneity and stratification of *NPM1*-mutated AML samples.


Methods: Using the MyAML® next generation sequencing (NGS) panel targeting 194 genes, we analyzed the genetic profile of 22 AML samples with driver mutations in the *NPM1* gene. The variants were annotated and interpreted following the Recommendation of the Association for Molecular Pathology, American Society of Clinical Oncology, College of American Pathologists, and American College of Medical Genetics and Genomics.

Results: In all cases the *NPM1* mutations co-occurred with mutations in the *FLT3* tyrosine kinase, *DNMT3A* methyltransferase, *IDH1* or *IDH2* isocitrate dehydrogenase genes. All *FLT3* mutations were located in the protein kinase catalytic domain and most mutations were internal tandem duplications (ITD). Four samples included more than one *FLT3*-ITD. Most of the *DNMT3A* mutations were detected at the mutational hotspot Arg882; all *IDH1* mutations were at the mutational hotspot Arg132; most *IDH2* mutations were at mutational hotspot Arg140, and all *SRSF2* mutations were found at mutational hotspot Pro95. Mutations in *IDH2* were more frequent than in *IDH1*, and mutations in these two genes did not co-occur in the same sample. *FLT3*-ITDs had the largest variability in allele frequency, when compared to mutations in *DNMT3A* and *IDH1/2*. The *NPM1*^{mut}/*FLT3*^{ITD}/*DNMT3A*^{mut} genotype, associated with poor prognosis in AML patients, was observed in 6 samples (~27%), close to a previously reported frequency³. All samples included mutations in the epigenetic modifiers *DNMT3A*, *IDH1/2*, *TET2* and *WT1*, which are involved in the DNA methylation/hydroxymethylation pathway⁶.

Mutational Profiles of AML Samples



Variant Allele Frequencies



ELN Recommendations for AML Treatment Based on Genetics

European LeukemiaNet - Conventional Care Regimens ⁵	
Consolidation therapy^{‡,§}	
<i>Younger patients (18-60/65 y)</i>	
Favorable-risk genetics	2-4 cycles of IDAC (1000-1500 mg/m ² IV over 3 h q12h, d1-3; or 1000-1500 mg/m ² IV over 3 h d1-5 or 6)
Intermediate-risk genetics	Allogeneic HCT from matched-related or unrelated donor
Adverse-risk genetics	2-4 cycles of IDAC (1000-1500 mg/m ² IV over 3 h q12h, d1-3; or 1000-1500 mg/m ² IV over 3 h d1-5 or 6), or High-dose therapy and autologous HCT
<i>Older patients (> 60/65 y)</i>	
Favorable-risk genetics	2-3 cycles of IDAC (500-1000 mg/m ² IV over 3 h q12h, d1-3; or 500-1000 mg/m ² IV over 3 h d1-5 or 6)
Intermediate/adverse-risk genetics	No established value of intensive consolidation therapy; consider allogeneic HCT in patients with low HCT-Comorbidity Index, or investigational therapy
Patients considered not candidates for intensive chemotherapy	
Low-dose cytarabine**	Low-dose cytarabine (20 mg q12h, SC, d1-10, q4 wk; until progression); not recommended in patients with adverse-risk genetics

‡ Patients, at least those aged 18 to 60 y, with newly diagnosed AML and activating *FLT3* mutations may be considered to receive additional therapy with midostaurin (administered after the chemotherapy). § Results from assessment of MRD should be taken into account for selecting consolidation therapy. ** In some countries used in a dosage of 20 mg/m² SC once daily. HCT, hematopoietic cell transplantation

Effect of Mutations in Mediators of the DNA Methylation Pathway




IDH mutations lead to production of 2-hydroxyglutarate (2HG), which inhibits TET2. TET2 mutations disrupt interaction with WT1. WT1 mutations cause failure to properly direct TET2 to target sites, either by disruption of the interaction itself or of DNA binding. aKG, alpha ketoglutarate⁶.

Conclusions

- The MyAML® NGS test allows a comprehensive, rapid, and cost-effective molecular profiling of AML samples.
- We confirmed the previously reported occurrences and statistically significant associations of driver mutations in AML samples^{2,3,8,11}. *NPM1* mutations associate with mutations in the DNA methylation/hydroxymethylation pathway and *FLT3*. Aberrant DNA methylation is a crucial process during tumorigenesis.
- The complexity of molecular interactions and their role in AML is a major challenge for clinical applications of genetic data. The classification "AML with *NPM1* mutation" only partially contributes to the risk stratification, based on the 2017 Recommendations by NCCN⁴ and ELN⁵, as well as on literature analysis. This classification does not inform about currently available treatments and clinical trials, which depend on the mutational profile of the AML sample obtained from sequencing of multiple genes.
- Therapies that combine conventional treatment and drugs targeting driver mutations are a promising strategy for improving the disease outcome. A combination of DNA methyltransferase and *FLT3* inhibitors have already been shown to be effective in AML patients⁹.
- NPM1*-positive minimal residual disease (MRD) has recently emerged as a sole prognostic factor for relapse, regardless of the presence of other gene alterations⁸. MRD monitoring allows precise and sensitive tracking of leukemia clones in order to apply personalized medicine⁷. Patients whose leukemia evolves from a pre-leukemic clone (e.g. *DNMT3A*, *TET2* mutations)¹⁰ may require therapy that exceeds standard chemotherapy. Monitoring of the pre-leukemic clone for acquisition of driver mutations (e.g. *NPM1*, *FLT3*^{ITD}) during remission might prompt initiation of aggressive therapy, such as allogeneic hematopoietic stem cell transplantation¹².

Combinational Therapies: Commercial and Developmental Drugs



A. New agents targeting epigenetic processes.¹³
 B. Rydapt (midostaurin) is used for treatment of AML patients who are *FLT3* mutation-positive.

Risk Stratification of AML Samples Based on Genetics

No	Genes	Guidelines	Scientific Publications (examples)
	<i>NPM1</i> <i>FLT3</i> <i>DNMT3A</i> <i>TET2</i>	NCCN ⁴ ELN ⁵	<i>NPM1</i> / <i>FLT3</i> / <i>DNMT3A</i> ³ <i>DNMT3A</i> / <i>FLT3</i> ¹⁴ <i>DNMT3A</i> / <i>TET2</i> ¹⁴
1	A	favorable	worse OS
18	A	favorable	worse DFS
21	D	intermediate	worse OS
5	A	intermediate	worse OS
7	A	intermediate	worse OS
13	A	intermediate	worse DFS
3	A	favorable	better OS
8	A	favorable	better OS
11	A	favorable	better OS
14	A	favorable	better OS
16	A	favorable	better OS
22	I	intermediate	better OS
9	A	favorable	worse DFS
12	A	favorable	worse DFS
2	A	favorable	
15	A	favorable	
19	A	favorable	better DFS
20	A	favorable	better DFS
4	A	favorable	better DFS
10	A	favorable	better DFS

References

- Gerstung M et al. Precision oncology for acute myeloid leukemia using a knowledge bank approach. *Nat Genet.* 2017 Mar; 49(3):332-340. PMID: 28092685.
- Bullinger L et al. Genomics of Acute Myeloid Leukemia: diagnosis and pathways. *J Clin Oncol.* 2017 Mar 20;35(9):934-946. PMID: 28297624.
- Papaemmanuil E et al. Genomic classification and prognosis in acute myeloid leukemia. *N Engl J Med.* 2016 Jun 9;374(23):2209-21. PMID: 27276561.
- NCCN Clinical Practice Guidelines in AML, Version 1.2017.
- Döhner H et al. Diagnosis and management of AML in adults: 2017 ELN recommendations from an international expert panel. *Blood.* 2017 Jan 26;129(4):424-447. PMID: 27895058.
- Rampal R, Figueroa ME. Wilms tumor 1 mutations in the pathogenesis of acute myeloid leukemia. *Haematologica.* 2016 Jun;101(6):672-9. PMID: 27252512.
- Cruz NM et al. Minimal residual disease in acute myelogenous leukemia. *Int J Lab Hematol.* 2017 May;39 Suppl 1:53-60. PMID: 28447422.
- Ivey A et al. Assessment of Minimal Residual Disease in Standard-Risk AML. *N Engl J Med.* 2016 Feb 4;374(5):422-33. PMID: 26789727.
- Ravandi F et al. Phase 2 study of azacitidine plus sorafenib in patients with acute myeloid leukemia and *FLT3* internal tandem duplication mutation. *Blood.* 2013 Jun 6;121(23):4655-62. PMID: 23613521.
- Jaiswal S et al. Age-related clonal hematopoiesis associated with adverse outcomes. *N Engl J Med.* 2014 Dec 25;371(26):2488-98. PMID: 25426837.
- Klco JM et al. Association Between Mutation Clearance After Induction Therapy and Outcomes in Acute Myeloid Leukemia. *JAMA.* 2015 Aug 25;314(8):811-22. PMID: 26305651.
- Koeffler HP, Leong G. Preleukemia: one name, many meanings. *Leukemia.* 2017 Mar;31(3):534-542. PMID: 27899806.
- Montalban-Bravo G, Garcia-Manero G. Myelodysplastic Syndromes – The Epigenetic Model for Drug Development? *Eur Oncol & Haematol.* 2016;12(1):13-4.
- Wang M et al. Molecular Mutations and Their Cooccurrences in Cytogenetically Normal Acute Myeloid Leukemia. *Stem Cells Int.* 2017;2017:6962379. PMID: 28197208.